



Pluslife FCM Card

One Nucleic Acid Test Card
Three Feline Pathogen Identification

FHV-1

C.felis

M.felis



For Veterinary Use Only

Features




TAT within 35 mins
FAST

Hands-on time in 1~2 mins
EASY

Comparable to PCR
ACCURATE

Up to 3 pathogens in 1 test
POWERFUL

Clinical Proven



Sensitivity: FHV-1 **100%**, C,felis **100%**, M.felis **95%**
Specificity: FHV-1 **100%**, C,felis **100%**, M.felis **100%**




Pluslife FCM Card

One Nucleic Acid Test Card
Three Feline Pathogen Identification


Pluslife Mini Dock is a rapid molecular diagnostic system that can be used to detect Feline herpesvirus I (FHV-1), Chlamydia felis (C.felis) and Mycoplasma felis (M.felis) in one single test. Our self-developed isothermal amplification technology, RHAM, provides excellent results comparable to qPCR tests. The product is compact and portable, with easy-to-operate procedure to obtain results within 35 minutes.

Test Procedure


1. Collect Sample



2. Mix Sample



3. Treat Sample



4. Run Test




Result Output



✓ Mini Dock



✓ PC Software



✓ Mobile App

Application



Pet Hospital/Clinic




Pet Shop



Animal CDC

Specification



Pluslife Mini Dock

Dimension	101 x 91 x 65 mm
Weight	210 g
Connection	Mobile phone or computer app
Operate temperature	15°C–28°C
Operate humidity	≤70%

Pluslife FCM Card

Target pathogen	FHV-1, C.felis, M.felis
Turnaround time (TAT)	35 mins
Shelf life	13 months at room temperature
Sample type	Eye/nasal/ oropharyngeal swab

Guangzhou Pluslife Biotech Co.,Ltd

 Floor 3, Block E, Runhui Science & Technology Park, No. 18 Shenzhou RD, Huangpu District, Guangzhou, Guangdong, China

 +86-20-31703986  market@pluslife.com  www.pluslife.com

